

RECEIVED  
CENTRAL FAX CENTER

MAR 30 2007

**Amendments to the Claims:**

This listing of claims will replace all prior versions, and listings, of claims in the application:

We claim:

**Claim 1 (Presently Amended)** An article of manufacture for care of contact lenses, comprising:

a container formed from poly (ethylene terephthalate) containing a composition comprising one or more surfactants and an ~~Alexidine~~-antimicrobial agent comprising 1,1' hexamethylene-bis[5-(2-ethylhexyl)biguanide.

**Claim 2 (Presently Amended)** The ~~composition~~article of claim 1 wherein said one or more surfactants are poloxamer surfactants or poloxamine surfactants or combinations thereof.

**Claim 3 (Presently Amended)** The ~~composition~~article of claim 1, wherein the composition further comprises at least one member selected from the group consisting of a buffering agent, a chelating agent, an osmolarity adjusting agent, and a surfactant having a HLB of 18 or above.

**Claim 4 (Presently Amended)** The ~~composition~~article of claim 1, wherein said ~~Alexidine~~ antimicrobial agent is present in an amount effective to disinfect a contact lens.

**Claim 5 (Presently Amended)** The ~~composition~~article of claim 1 wherein the composition comprises about 0.1 to about 6.0 weight percent of said surfactant and about 0.05 to about 0.5 weight percent of said antimicrobial agent.

**Claim 6 (Presently Amended)** The ~~composition~~article of claim 1 wherein the composition further comprises:

a chelating agent and a buffering agent selected from the group consisting borate buffers, phosphate buffers and citrate buffers.

**Claim 7 (Presently Amended)** The ~~composition~~article of claim 6, wherein the composition comprises at least one member selected from the group consisting of poloxamer and poloxamine surfactants having HLB values of 18 or greater.

**Claim 8 (Withdrawn)** A method of enhancing biocidal efficacy of a lens care solution comprising:

packaging a lens care solution containing one or more surfactants and one or more antimicrobial agents in a container formed from poly(ethylene terephthalate).

**Claim 9 (Withdrawn)** A method of enhancing stability of a lens care solution comprising:

packaging a lens care solution containing one or more surfactants and one or more antimicrobial agents in a container formed from poly(ethylene terephthalate).

**Claim 10 (Withdrawn)** A method of increasing shelf-life of a lens care solution comprising:

packaging a lens care solution containing one or more surfactants and one or more antimicrobial agents in a container formed from poly(ethylene terephthalate).

**Claim 11 (Withdrawn)** The method of claim 8, 9 or 10 wherein the solution further comprises at least one member selected from the group consisting of a buffering agent, a chelating agent, an osmolarity adjusting agent, and a surfactant having a HLB value of 18 or greater.

**Claim 12 (Withdrawn)** The method of claim 8, 9 or 10 wherein the solution further comprises an antimicrobial agent in an amount effective to disinfect a contact lens.

**Claim 13 (Withdrawn)** The method of claim 8, 9 or 10 wherein the solution comprises about 0.05 to about 0.5 weight percent of said antimicrobial agent.

**Claim 14 (Withdrawn)** The method of claim 8, 9 or 10 wherein the solution further comprises: a chelating agent and a buffering agent selected from the group consisting borate buffers, phosphate buffers and citrate buffers.

**Claim 15 (Withdrawn)** The method of claim 8, 9 or 10 wherein the composition further comprises a surfactant having a HLB value of 18 or greater.

**Claim 16 (Withdrawn)** The method of claim 8, 9 or 10 wherein the solution comprises at least one member selected from the group consisting of poloxamer and poloxamine surfactants having a HLB value of 18 or greater.

**Claim 17 (Withdrawn)** The method of claim 8, 9 or 10 wherein said one or more surfactants are selected from the group consisting of Pluronic P123<sup>TM</sup>, Pluronic L42<sup>TM</sup>, Pluronic L62<sup>TM</sup>, Pluronic L72<sup>TM</sup>, Pluronic L92<sup>TM</sup>, Pluronic P103<sup>TM</sup>, Pluronic R 12R3<sup>TM</sup>, Pluronic R 17R1<sup>TM</sup>, Pluronic R 17R2<sup>TM</sup>, Pluronic R 31R1<sup>TM</sup>, Pluronic R 31R2<sup>TM</sup>, Pluronic R 31R4<sup>TM</sup>, Tetronic 701<sup>TM</sup>, Tetronic 702<sup>TM</sup>, Tetronic 901<sup>TM</sup>, Tetronic 1101<sup>TM</sup>, Tetronic 1102<sup>TM</sup>, Tetronic 1301<sup>TM</sup>, Tetronic 1302<sup>TM</sup>, Tetronic 1501<sup>TM</sup>, Tetronic 1502<sup>TM</sup>, Tetronic R 50R1<sup>TM</sup>, Tetronic R 50R4<sup>TM</sup>, Tetronic R 70R1<sup>TM</sup>, Tetronic R 70R2<sup>TM</sup>, Tetronic R 70R4<sup>TM</sup>, Tetronic R 90R1<sup>TM</sup>, Tetronic R 90R4<sup>TM</sup>, Tetronic R 110R1<sup>TM</sup>, Tetronic R 110R2<sup>TM</sup>, Tetronic R 110R7<sup>TM</sup>, Tetronic R 130R1<sup>TM</sup>, Tetronic R 130R2<sup>TM</sup>, Tetronic R 150R1<sup>TM</sup>, Tetronic R 150R4<sup>TM</sup> and Tetronic R 150R8<sup>TM</sup>.